


Patient Biometrics System (PBS) Processor User Guide

Changelog Version 2301.02

1. Fix the bug: *Data source rejected establishment of connection, message from server: "Too many connections"*.
2. Enhance fingerprints processing feedback.
3. Enhanced export performance.
4. Fix application freeze on a single thread.
5. Flexibility to configure database parameters.
6. Launch this documentation from the user interface by clicking  button.

How-To

Step 1: Edit the **app.properties** with the parameters to match your NigeriaMRS host server.





Figure 1: app.properties configuration file

Step 2: Launch the **PBSprocessor.jar** application

Name	Type
docs	File folder
lib	File folder
app	PROPERTIES File
<u>PBSprocessor</u>	Executable Jar File

Figure 2 : PBS processor application launcher

Step 3: Click the  button to browse the exported PBS file and start processing the fingerprints.

Step 4: Browse to the folder location where the .json PBS exported file(s) is/are and click **Select Folder** from the folder dialog box as shown in figure 3.




Figure 3: PBS files browser




Figure 4: Fingerprints processed